

MicroClime™ Environmental Lid

**Features**

This lid is designed to fit SBS conforming multi-well plates, 127.76 mm x 85.48 mm. The lid carries stabilized liquid appropriate to the plates being covered. The lid dramatically reduces edge effects in wells including hydration and evaporation of solvents. Lids can be stacked on each other. Reservoir can be filled manually or robotically from long or short edges. Reservoir fluid does not contact wells.



**Specifications**

- Item #:** LL-0300
- Composition:** COC
- Felt material:** Polypropylene
- Format:** Fit standard SBS plates
- Color:** Clear
- Sterile:** Ask
- Reservoir volume:** 9mL
- Qty/Case:** 50 lids/case
- Qty/Package:** 10 lids/package

