

Aqueous Applications

384

Wells

Plate Type

LDV

(Low Dead Volume)

PP

(Polypropylene)

Fluid Type

Salt-based buffer
± nucleic acid,
cell culture media
w/o serum

Salt-based buffer
+ protein
Cell culture media
+ serum

Salt-based buffer
± protein, ± nucleic
acid; Cell culture
media ± serum

Surfactant-
containing
solutions

Volume Range (μL)

3-12

6-14

20-50

20-50

Working Range (μL)

9

8

30

30

Heat Seal
Compatible

YES

YES

YES

YES